




ANAESTHESIA

CATHETER FORCEPS



301.15
150mm / 6"
For Babies

302.20
200mm / 8"
For Children

303.25
250mm / 10"
For Adults

MAGILL
301.15 - 303.25

LARYNGOSCOPES



MILLER
 350.00 - 352.00



353.00



353.01



353.02



353.03



353.04

350.00 Laryngoscope Set consisting of:
 392.00 Batterie Handle (without batterie)
 353.01 Blade, Fig.1
 353.02 Blade, Fig.2
 353.03 Blade, Fig.3

351.00
 as 350.00 with 4 blades

352.00 Miller - Baby
 Laryngoscope Set for Baby and Infant,
 consisting of:
 392.00 Batterie Handle (without batterie)
 353.00 Blade, Fig.0
 353.01 Blade, Fig.1

Miller Laryngoscope Blades
 353.00 Fig.0, working lenght 55mm, for baby's
 353.01 Fig.1, working lenght 80mm, for children's
 353.02 Fig.2, working lenght 130mm, for adolescent's
 353.03 Fig.3, working lenght 170mm, for women's
 353.04 Fig.4, working lenght 180mm, for men's



392.00



395.00



DD 395.01